



## Product Data Sheet - 2 Liter Vacuum Filter Funnel System Part Number VSS – 1020

**Velocity Scientific Solutions**  
PO Box 1547, Canning Vale DC WA  
6970  
T: 1300 855 315  
F: 1300 855 316  
E: [info@velocityscientific.com.au](mailto:info@velocityscientific.com.au)  
W: [www.velocityscientific.com.au](http://www.velocityscientific.com.au)

### Description

This Vacuum Filter Funnel is made of durable polypropylene and stainless steel. It is easy to set up and use and is virtually indestructible. The unit consists of:

Qty.	Description
1 ea	2 Liter Funnel
1 ea	Filter Retainer Ring
1 ea	Stainless Steel Membrane Support
1 ea	Base Unit with Vacuum Port
1 ea	8-piece Adapter Kit


### Assembly

Assemble the unit for use as follows:

- Secure Base Unit (E) to receiving bottle. The receiving bottle must have a GL-45 threaded neck. The receiving bottle should never exceed 2 liters in volume. It is recommended that safety coated receiving bottles be used.
- Install the Stainless Steel Membrane Support (D) into the Base Unit (E). This will fit snugly.
- Apply the 47mm Membrane Filter (C) to rest on the Stainless Steel Membrane Support.
- Install the Filter Retainer Ring (B). This ring protects the Membrane Filter and enhances the seal between the Funnel and Base Unit. The unit may be used without the Filter Retainer Ring, but great care must be taken to avoid damaging the Membrane Filter.
- Seat the Funnel (A) on top of the Filter Retainer Ring (B) and raise the white cap to engage the threads on the Funnel. It is recommended that the cap be turned to tighten and seal the Funnel to the Base Unit, rather than turning the Funnel.
- Secure the unit with a ring stand. A three-finger clamp can be used to hold the base unit, securely fixing the entire filtration apparatus to the ring stand.
- Thread the 1/4"-28 to female Luer adapter into the vacuum port. Apply the appropriate barb to interface the unit to the vacuum source.
- Apply vacuum and use the filter.

### Safety Precautions

- Ensure fluid compatibility with the filter system materials.
- Always wear proper personal protective gear appropriate to the task.
- Always use a safety coated reservoir bottle in conjunction with the filter system.
- Never pull vacuum on a reservoir bottle larger than 2 liters.



**Figure 1** – Components of the Vacuum Filter Funnel and their proper orientation and assembly. **A** – Funnel, **B** – Filter Retainer Ring, **C** – 47mm Filter Membrane, **D** – Stainless Steel Membrane Support, **E** – Base Unit with Vacuum Port, **F** – GL-45 Cap Closure