

Environmental Vials, Closures and Septum



- Available in clear or amber borosilicate glass
- Choice of closure and liner, including patented Top Hat™ liner system
- Fixed (bonded) closures eliminate liner fallout

Vials - Class 1

Cat. No.	Description	Qty
D0376-20	20mL Clear Vial, 28x56mm, 24-400mm Thread	100 Pack
D0390-20	20mL Amber Vial, 28x56mm, 24-400mm Thread	100 Pack
320024-2856	20mL Clear Vial, 28x56mm, 24-400mm Thread	144 Pack
320024-2856A	20mL Amber Vial, 28x56mm, 24-400mm Thread	144 Pack
D0376-40	40mL Clear Vial, 28x95mm, 24-400mm Thread	100 Pack
D0390-40	40mL Amber Vial, 28x95mm, 24-400mm Thread	100 Pack
340024-2895	40mL Clear Vial, 28x95mm, 24-400mm Thread	144 Pack
340024-2895A	40mL Amber Vial, 28x95mm, 24-400mm Thread	144 Pack
D0376-60	60mL Clear Vial, 27x140mm, 24-400mm Thread	100 Pack
D0390-60	60mL Amber Vial, 27x140mm, 24-400mm Thread	100 Pack
36024-27140	60mL Clear Vial, 27x140mm, 24-400mm Thread	144 Pack
36024-27140A	60mL Amber Vial, 27x140mm, 24-400mm Thread	144 Pack

Closures - Class 1

Cat. No.	Description	Qty
5310-24	24-400mm Black, Polypropylene Open Hole Closure, 100 per pack	1,000 Case
5310-24W	24-400mm White, Polypropylene Open Hole Closure, 100 per pack	1,000 Case
5310-24414W	24-414mm White, Polypropylene Open Hole Closure, 100 per pack	1,000 Case
5320-24	24-400mm Solid Top, Black Polypropylene Unlined Closure, 100 per pack	1,000 Case
5320-24W	24-400mm Solid Top, White Polypropylene Unlined Closure, 100 per pack	1,000 Case
5340-24	24-400mm Black, Polypropylene Glass Filled Open Hole Closure, 100 per pack	1,000 Case

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Closures with Liners - Class 1

Cat. No.	Description	Qty
34-527	24-414mm, White Polypropylene Open Hole Closure, Fixed 0.125" PTFE/Silicone Liner	100 Pack
34-527/144	24-414mm, White Polypropylene Open Hole Closure, Fixed 0.125" PTFE/Silicone Liner	144 Pack
34-527/M	24-414mm, White Polypropylene Open Hole Closure, Fixed 0.125" PTFE/Silicone Liner	1,000 Pack
34-375BLK	24-414mm, Black Polypropylene Open Hole Closure, Fixed 0.060" PTFE/Silicone Liner	100 Pack
D0396-24	24-400mm Black Phenolic Solid Top Closure, PE Cone	100 Pack
D0399-24	24-400mm, White Polypropylene Solid Top Closure, PTFE Lined	100 Pack
D0399-24/144	24-400mm, White Polypropylene Solid Top Closure, PTFE Lined	144 Pack
D0400-24	24-400mm, White Polypropylene Solid Top Closure, F217 Lined	100 Pack
D0400-24/144	24-400mm, White Polypropylene Solid Top Closure, F217 Lined	144 Pack
806550TH-24	24-400mm Black Polypropylene Open Top Closure, .0125" Top Hat™ PTFE/Silicone Lined (Patented)	100 Pack

Septum - Class 1

Cat. No.	Description	Qty
38-140	24mm x 0.015" PTFE	100 Pack
610050-24	24mm x 0.100" PTFE/Silicone Septa, 100 per pack	1,000 Case
612550-24	24mm x 0.125" PTFE/Silicone Septa, 100 per pack	1,000 Case

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*