



Biotage Cross Reference

Santai Product	Santai Item#	Description	Biotage Product	Type of Column	Part Number
Standard Series (Surface area 500m ² /g, pH 6.5-7.5) 	S-5101-0004	4g flash cartridge, Irregular silica gel 40-63 μm 60A	Biotage ZIP Cartridge	Biotage ZIP cartridge, 5 g	440-0500-DZ-20
	S-5101-0012	12g flash cartridge, Irregular silica gel 40-63 μm 60A		Biotage ZIP cartridge, 10 g	440-1000-EZ-20
	S-5101-0025	25g flash cartridge, Irregular silica gel 40-63 μm 60A		Biotage ZIP cartridge, 30 g	440-3000-FZ-20
	S-5101-0040	40g flash cartridge, Irregular silica gel 40-63 μm 60A		Biotage ZIP cartridge, 45 g	440-4500-SZ-20
	S-5101-0080	80g flash cartridge, Irregular silica gel 40-63 μm 60A		Biotage ZIP cartridge, 80 g	440-8000-JZ-20
	S-5101-0120	120g flash cartridge, Irregular silica gel 40-63 μm 60A		Biotage ZIP cartridge, 120 g	440-120G-UZ-20
	S-5101-0220	220g flash cartridge, Irregular silica gel 40-63 μm 60A		/	/
	S-5101-0330	330g flash cartridge, Irregular silica gel 40-63 μm 60A		/	/
	S-5101-0800	800g flash cartridge, Irregular silica gel 40-63 μm 60A		/	/
	S-5101-1600	1600g flash cartridge, Irregular silica gel 40-63 μm 60A		/	/
	S-5101-3000	3000g flash cartridge, Irregular silica gel 40-63 μm 60A		/	/

Santai Product	Santai Item#	Description	Biotage Product	Type of Column	Part Number
Standard Series (Surface area 680m ² /g, pH 6.0-8.0) 	S-2101-0004-SP(H)	4g flash cartridge, High-capacity Spherical silica gel 50μm 54A	Biotage ZIP Spherical Cartridge	ZIP Spherical Silica cartridge, 5 g	445-0500-DZ-20
	S-2101-0012-SP(H)	12g flash cartridge, High-capacity Spherical silica gel 50μm 54A		ZIP Spherical Silica cartridge, 10 g	445-1000-EZ-20
	S-2101-0025-SP(H)	25g flash cartridge, High-capacity Spherical silica gel 50μm 54A		ZIP Spherical Silica cartridge, 30 g	445-3000-FZ-20
	S-2101-0040-SP(H)	40g flash cartridge, High-capacity Spherical silica gel 50μm 54A		ZIP Spherical Silica cartridge, 45 g	445-4500-SZ-20
	S-2101-0080-SP(H)	80g flash cartridge, High-capacity Spherical silica gel 50μm 54A		ZIP Spherical Silica cartridge, 80 g	445-8000-JZ-20
	S-2101-0120-SP(H)	120g flash cartridge, High-capacity Spherical silica gel 50μm 54A		ZIP Spherical Silica cartridge, 120 g	445-120G-UZ-20
	S-2101-0220-SP(H)	220g flash cartridge, High-capacity Spherical silica gel 50μm 54A		/	/
	S-2101-0330-SP(H)	330g flash cartridge, High-capacity Spherical silica gel 50μm 54A		/	/
	S-2101-0800-SP(H)	800g flash cartridge, High-capacity Spherical silica gel 50μm 54A		/	/
	S-2101-1600-SP(H)	1600g flash cartridge, High-capacity Spherical silica gel 50μm 54A		/	/
	S-2101-3000-SP(H)	3000g flash cartridge, High-capacity Spherical silica gel 50μm 54A		/	/