

2.0mL Wide Mouth 10-425 Screw Vials, Inserts and Caps



Vial Volume: 2.0 mL
Width: 12mm
Height: 32mm
Opening Diameter: 6mm
Cap Style: 10/425 Screw Cap

40% larger opening prevents broken needles due to increased target area. Choose from clear or amber Type I borosilicate glass. ID design incorporates the StepVial feature that precisely centers a limited volume insert in the vial neck. Available with graduated marking spots which can be written on with a #2 pencil.

Agilent	Perkin Elmer	Varian	Waters	Velocity	Description
5183-4401	392611646		WAT063300	32010-1232	12x32mm Clear Wide Mouth Vial, 10mm/425 Thread
			WAT270946	806550T-1232	12x32mm Clear Vial with Red PTFE/Silicone Septum and Cap
	392611650		WAT094174	806050-10	10/425 Black PP Open-Hole Cap with PTFE/Silicone Septum
			WAT058876	806070-10	10/425 Black PP Open-Hole Cap with PTFE/Silicone Septa with Slit
			WAT058875	5310-10	10/425 Black Polypropylene (PP) Open-Hole Cap
			WAT058874	606050-10	10mm x 0.060" Red PTFE/Silicone Septum
5183-2085		392611596		4025-629	250uL LV Glass Insert for Wide Mouth Vial
5181-1270		392611594	WAT094170	4025BS-629	250uL LV Glass Insert with PE Bottom Spring for Wide Mouth Vial
5181-3377				4035-630	350uL LV Glass Flat Bottom Insert for Wide Mouth Vial
5182-0549				4025PBS-631	250uL LV PP Insert with PE Bottom Spring for Wide Mouth Vial
5183-2087				4035P-630	350uL LV PP Flat Bottom Insert for Wide Mouth Vial



30210-1232

Vials	
Description	Part Number
10-425 Standard Mouth Screw Top Vials	
2.0mL Clear 10-425 Wide Mouth Vial 100 Pack	32010-1232
2.0mL Clear 10-425 Wide Mouth Vial 100 Pack - Deactivated	32010-1232Z
2.0mL Amber 10-425 Wide Mouth Vial 100 Pack	32010-1232A
2.0mL Amber 10-425 Wide Mouth Vial 100 Pack - Deactivated	32010-1232AZ
2.0mL Clear 10-425 Wide Mouth Vial with White Graduation Spot 100 Pack	32010E-1232
2.0mL Amber 10-425 Wide Mouth Vial with White Graduation Spot 100 Pack	32010E-1232A
10-425 Limited Volume Vials	
300µL Glass Insert in clear Glass 10-425 Interlocked™ Vial 100 Pack	30210-1232
300µL Glass Insert in Amber 10-425 Glass Interlocked™ Vial 100 Pack	30210-1232A
300µL Glass Insert in Clear 10-425 Glass Interlocked™ Vial with Marking Spot 100 Pack	30210M-1232
300µL Glass Insert in Amber 10-425 Interlocked™ Vial with Marking Spot 100 Pack	30210M-1232A
Plastic 10-425 Standard Mouth Screw Top Vials	
500µL 10-425 Wide Mouth Polypropylene Vial 100 Pack	30510P-1232
500µL Clear 10-425 Wide Mouth Polypropylene Vial 100 Pack	30510CP-1232
500µL TPX 10-425 Wide Mouth (Polymethylpentene) Vial 100 Pack	30510T-1232
750µL 10-425 Wide Mouth Polypropylene Vial 100 Pack	30710P-1232
750µL Clear 10-425 Wide Mouth Polypropylene Vial 100 Pack	30710CP-1232
750µL TPX 10-425 Wide Mouth (Polymethylpentene) Vial 100 Pack	30710T-1232



30510P-1232, 30510CP-1232, 30510T-1232



30710P-1232, 30710CP-1232, 30710T-1232

2.0mL Wide Mouth 10-425 Screw Vials, Inserts and Caps



Open Hole Caps (No Liner)

Description	Part Number
Black, 10-425, Polypropylene Large Open Hole Cap 100 Pack	5310-10
White, 10-425, Polypropylene Large Open Hole Cap 100 Pack	5310-10W
Yellow, 10-425, Polypropylene Large Open Hole Cap 100 Pack	5310-10Y
Blue, 10-425, Polypropylene Large Open Hole Cap 100 Pack	5310-10B
Red, 10-425, Polypropylene Large Open Hole Cap 100 Pack	5310-10R




Solid Top Storage Caps

Description	Part Number
Solid Top, Black Polypropylene Unlined 100 Pack	5320-10
Solid Top, Black Polypropylene, PTFE/F217 Lined 100 Pack	5360-10



10mm Septa/Liners

Description	Part Number
10mm x 0.010" White PTFE Septa 100 Pack	601010-10
10mm x 0.040" Red PTFE/Silicone/Red PTFE Septa 100 Pack	604060-10
10mm x 0.040" Red PTFE/Silicone/Red PTFE w/Slit Septa 100 Pack	604075-10
10mm x 0.045" PTFE/Red Rubber Septa 100 Pack	604040-10
10mm x 0.060" Red PTFE/Silicone w/Slit Septa 100 Pack	606070-10
10mm x 0.065" Red PTFE/Silicone Septa 100 Pack	606550-10



Preassembled 10-425 Caps and Septa

Description	Part Number
Black Closure, PTFE/Red Rubber Septa, 0.040" 100 Pack	804040-10
Black Closure, PTFE/Silicone/PTFE Septa, 0.040" 100 Pack	804060-10
Black Closure, Red PTFE/Silicone Septa, 0.060" w/Slit 100 Pack	806070-10
Black Closure, Red PTFE/Silicone Septa, 0.065" 100 Pack	806550-10



4025BS-629, 401TS-529

Inserts

Description	Part Number
Conical Glass Inserts	
50µL Glass Insert with Bottom Spring 100 Pack	4005BS-625
50µL Glass Insert with Bottom Spring 100 Pack - Deactivated	4005BS-625Z
100uL Glass Insert with Glass Flange 100 Pack	401TF-529
100uL Glass Insert with Top Spring 100 Pack	401TS-529
250µL Glass Big Mouth Insert, No Spring 100 Pack	4025-629
250µL Glass Big Mouth Insert, No Spring 100 Pack - Deactivated	4025-629Z
250µL Glass Big Mouth Insert with Bottom Spring 100 Pack	4025BS-629
250µL Glass Big Mouth Insert with Bottom Spring 100 Pack - Deactivated	4025BS-629Z
Flat Bottom Glass Inserts	
350µL Glass Flat Bottom Insert 100 Pack	4035-630
350µL Glass Flat Bottom Insert 100 Pack - Deactivated	4035-630Z
350µL Glass Flat Bottom Insert with ID Ring 100 Pack	4035-630ID
350µL Glass Flat Bottom Insert with ID Ring 100 Pack - Deactivated	4035-630IDZ
Conical Polypropylene Inserts	
100µL Polypropylene Insert with Top Spring 100 Pack	401PTS-530
250µL Polypropylene Big Mouth Insert, No Spring 100 Pack	4025P-631
250µL Polypropylene Big Mouth Insert with Bottom Spring 100 Pack	4025PBS-631



4035-630ID, 4035-630