

2.0mL Standard Mouth Crimp Top Vials, Inserts and Caps



Vial Volume: 2.0 mL
Width: 12mm
Height: 32mm
Opening Diameter: 4mm
Cap Style: 11mm Crimp

Neck finish allows use of Snap Top Caps™, patented Poly Crimp™ Seals or standard aluminum seals. Snap Ring™ finish eliminates the need for crimping or decapping. Designed to work in Agilent and other autosamplers. Choose from clear or amber Type I borosilicate glass. ID design incorporates the StepVial feature that precisely centers a limited volume insert in the vial neck. Available with graduated marking spots which can be written on with a #2 pencil.

Agilent	Perkin Elmer	Varian	Velocity	Description
	N930-1385	392611626	32011-1232	12x32mm Clear Standard Mouth Vial, 11mm Crimp
	N930-2680	392611628	32011-1232A	12x32mm Amber Standard Mouth Vial, 11mm Crimp



30211M-1232, 30211-1232

Vials	
Description	Part Number
Standard Vials	
2.0mL Clear Standard Crimp Top Vial 100 Pack	32011-1232
2.0mL Clear Standard Crimp Top Vial 100 Pack - Deactivated	32011-1232Z
2.0mL Amber Standard Crimp Top Vial 100 Pack	32011-1232A
2.0mL Amber Standard Crimp Top Vial 100 Pack - Deactivated	32011-1232AZ
2.0mL Clear Standard Crimp Top Vial with White Graduation Spot 100 Pack	32011E-1232
2.0mL Amber Standard Crimp Top Vial with White Graduation Spot 100 Pack	32011E-1232A
Interlocked Vials with Fused glass Inserts	
300µL Clear Interlocked™ Vial/Insert 100 Pack	30211-1232
300µL Amber Interlocked™ Vial/Insert 100 Pack	30211-1232A
300µL Clear Interlocked™ Vial/Insert with Marking Spot 100 Pack	30211M-1232
300µL Amber Interlocked™ Vial/Insert with Marking Spot 100 Pack	30211M-1232A



11mm Aluminium Caps

Description	Part Number				
	Silver	Blue	Gold	Red	Green
11mm Seal, 10mil PTFE 100 Pack	5110-11	-	-	-	-
11mm Seal, Clear PTFE/Red Rubber 100 Pack	5140-11	5140-11B	5140-11Y	5140-11R	5140-11G
11mm Seal, PTFE/Silicone 100 Pack	5150-11	5150-11B	5150-11Y	5150-11R	5150-11G
11mm Seal, PTFE/Silicone/PTFE 100 Pack	5160-11	5160-11B	5160-11Y	5160-11R	5160-11G
11mm Seal, PTFE/Silicone w/Slit 100 Pack	5170-11	-	-	-	-



4005BS-425, 401BS-530, 401F-529

Inserts	
Description	Part Number
Conical Inserts	
50µL Conical Glass Insert w/Bottom Spring 100 Pack	4005BS-425
50µL Conical Glass Insert w/Bottom Spring 100 Pack - Deactivated	4005BS-425Z
100µL Glass Insert w/Plastic Flange 100 Pack	401F-529
100µL Glass Insert w/Bottom Spring 100 Pack	401BS-530
100µL Glass Insert w/Bottom Spring 100 Pack - Deactivated	401BS-530Z
100µL Glass Insert, No Spring 100 Pack	401-531
100µL Glass Insert, No Spring 100 Pack - Deactivated	401-531Z
150µL Glass Insert, w/Bottom Spring 100 Pack	4015BS-529
150µL Glass Insert, w/Bottom Spring 100 Pack - Deactivated	4015BS-529Z
100µL Polypropylene Insert, No Spring 100 Pack	401P-530
100µL Polypropylene Insert w/Bottom Spring 100 Pack	401PBS-530
Flat Bottom Inserts	
250µL Glass Flat Bottom Insert 100 Pack	4025-531
250µL Glass Flat Bottom Insert 100 Pack - Deactivated	4025-531Z
250µL Glass Flat Bottom Insert w/ID Ring 100 Pack	4025-531ID
250µL Glass Flat Bottom Insert w/ID Ring 100 Pack - Deactivated	4025-531IDZ
250µL Polypropylene Flat Bottom Insert 100 Pack	4025P-531



4025-531ID, 4025-531