

2.0mL Wide Mouth Crimp Top Vials, Inserts and Caps



Vial Volume: 2.0 mL
Width: 12mm
Height: 32mm
Opening Diameter: 6mm
Cap Style: 11mm Crimp

40% larger opening prevents broken needles due to increased target area. Neck finish allows use of patented Poly Crimp™ Seals or standard aluminium seals. Designed to work in Agilent and other autosamplers. Choose from clear or amber Type I borosilicate glass. ID design incorporates the StepVial feature that precisely centres a limited volume insert in the vial neck. Available with graduated marking spots which can be written on with a #2 pencil.

Agilent	Varian	Waters	Velocity	Description
5181-3375	392611634	WAT094222	32011L-1232	12x32mm Clear Wide Mouth Vial, 11mm Crimp
5182-0543			32011LE-1232	12x32mm Clear Wide Mouth Vial with Marking Spot, 11mm Crimp
5181-3376	392611635	WAT094223	3211LE-1232A	12x32mm Amber Wide Mouth Vial with Marking Spot, 11mm Crimp
5182-3454	186000984		31511-1232	12x32mm High Recovery Vial (HRV), 1.5mL, Wide Mouth Crimp
	392611875		5540-11	11mm Poly Crimp Seal with PTFE/Butyl Septum
	392611873		5550-11	11mm Poly Crimp Seal with PTFE/Silicone Septum
	392611874		5570-11	11mm Poly Crimp Seal with PTFE/Silicone Septum with Slit
9301-0977			30111G-1232	100µL Glass Insert in Clear Polypropylene Vial
5183-2085	392611596		4025-629	250uL LV Glass Insert for Wide Mouth Vial
5181-1270	392611594	WAT094170	4025BS-629	250uL LV Glass Insert with PE Bottom Spring for Wide Mouth Vial
5181-3377			4035-630	350uL LV Glass Flat Bottom Insert for Wide Mouth Vial
5182-0549			4025PBS-631	250uL LV PP Insert with PE Bottom Spring for Wide Mouth Vial
5183-2087			4035P-630	350uL LV PP Flat Bottom Insert for Wide Mouth Vial



31511-1232

Vials	Description	Part Number
11mm Big Mouth Crimp Vials		
	1.5mL Clear Big Mouth High Recovery Vial 100 Pack	31511-1232
	2.0mL Clear Big Mouth Crimp Top Vial 100 Pack	32011L-1232
	2.0mL Clear Big Mouth Crimp Top Vial 100 Pack - Deactivated	32011L-1232Z
	2.0mL Amber Big Mouth Crimp Top Vial 100 Pack	32011L-1232A
	2.0mL Amber Big Mouth Crimp Top Vial 100 Pack - Deactivated	3211L-1232AZ
	2.0mL Clear Big Mouth Crimp Top Vial with White Graduation Spot 100 Pack	32011LE-1232
	2.0mL Amber Big Mouth Crimp Top Vial with White Graduation Spot 100 Pack	32011LE-1232A
11mm Big Mouth Limited Volume Crimp Vials		
	100µL Glass Insert in Clear Polypropylene Vial 100 Pack	30111G-1232
	100µL Glass Insert in Amber Polypropylene Vial 100 Pack	30111G-1232A
	300µL Clear Interlocked™ Vial/Insert 100 Pack	30211L-1232
	300µL Amber Interlocked™ Vial/Insert 100 Pack	30211L-1232A
	300µL Clear Interlocked™ Vial/Insert with Marking Spot 100 Pack	30211LM-1232
	300µL Amber Interlocked™ Vial/Insert with Marking Spot 100 Pack	30211LM-1232A
	600µL Glass Insert in Clear Polypropylene Vial 100 Pack	30611G-1232
Plastic Crimp/Snap Top Vials		
	100µL Clear Polypropylene Crimp/Snap Vial 100 Pack	30111CP-1232
	500µL Clear Polypropylene Crimp/Snap Vial 100 Pack	30511CP-1232
	750µL Clear Polypropylene Crimp/Snap Vial 100 Pack	30711CP-1232



30211LM-1232, 30211L-1232



30611G-1232



30511CP-1232, 30111CP-1232

2.0mL Wide Mouth Crimp Top Vials, Inserts and Caps



Aluminium Seals (11mm) - Requires an 11mm Crimper					
Description	Part Number				
	Silver	Blue	Yellow	Red	Green
11mm Seal, 10mil PTFE 100 Pack	5110-11	-	-	-	-
11mm Seal, Clear PTFE/Red Rubber 100 Pack	5140-11	5140-11B	5140-11Y	5140-11R	5140-11G
11mm Seal, PTFE/Silicone 100 Pack	5150-11	5150-11B	5150-11Y	5150-11R	5150-11G
11mm Seal, PTFE/Silicone/PTFE 100 Pack	5160-11	5160-11B	5160-11Y	5160-11R	5160-11G
11mm Silver Seal, PTFE/Silicone w/Slit 100 Pack	5170-11	-	-	-	-



Poly Crimp Seals (11mm) - Metal Free Seal, Requires an 11mm Crimper					
Description	Part Number				
	Clear	Blue	Yellow	Red	Green
11mm Seal, 10mil PTFE 100 Pack	5510-11	5510-11B	5510-11Y	5510-11R	5510-11G
11mm Seal, Clear PTFE/Red Rubber 100 Pack	5540-11	5540-11B	5540-11Y	5540-11R	5540-11G
11mm Seal, PTFE/Silicone 100 Pack	5550-11	5550-11B	5550-11Y	5550-11R	5550-11G
11mm Seal, PTFE/Silicone/PTFE 100 Pack	5560-11	5560-11B	5560-11Y	5560-11R	5560-11G
11mm Silver Seal, PTFE/Silicone w/Slit 100 Pack	5570-11	5570-11B	5570-11Y	5570-11R	5570-11G



4005BS-625, 4025BS-529, 401TS-529

Inserts	
Description	Part Number
Conical Glass Inserts	
50µL Glass Insert with Bottom Spring 100 Pack	4005BS-625
50µL Glass Insert with Bottom Spring 100 Pack - Deactivated	4005BS-625Z
100uL Glass Insert with Glass Flange 100 Pack	401TF-529
100uL Glass Insert with Top Spring 100 Pack	401TS-529
250µL Glass Big Mouth Insert, No Spring 100 Pack	4025-629
250µL Glass Big Mouth Insert, No Spring 100 Pack - Deactivated	4025-629Z
250µL Glass Big Mouth Insert with Bottom Spring 100 Pack	4025BS-629
250µL Glass Big Mouth Insert with Bottom Spring 100 Pack - Deactivated	4025BS-629Z
Flat Bottom Glass Inserts	
350µL Glass Flat Bottom Insert 100 Pack	4035-630
350µL Glass Flat Bottom Insert 100 Pack - Deactivated	4035-630Z
350µL Glass Flat Bottom Insert with ID Ring 100 Pack	4035-630ID
350µL Glass Flat Bottom Insert with ID Ring 100 Pack - Deactivated	4035-630IDZ
Conical Polypropylene Inserts	
100µL Polypropylene Insert with Top Spring 100 Pack	401PTS-530
250µL Polypropylene Big Mouth Insert, No Spring 100 Pack	4025P-631
250µL Polypropylene Big Mouth Insert with Bottom Spring 100 Pack	4025PBS-631



4035-630ID, 4035-630



4025P-631, 4025PBS-631

Looking for Storage Racks and Crimpers?
See page 123 for more details



4025BS-629, 401TS-529