

## 2.0mL Wide Mouth Snap Top Vials, Inserts and Caps



**Vial Volume:** 2.0 mL  
**Width:** 12mm  
**Height:** 32mm  
**Opening Diameter:** 6mm  
**Cap Style:** 11mm Crimp/Snap

Neck finish allows use of Snap Top Caps™, patented Poly Crimp™ Seals or standard aluminium seals. Snap Ring™ finish eliminates the need for crimping or decapping. Designed to work in Agilent and other autosamplers. Choose from clear or amber Type I borosilicate glass. ID design incorporates the StepVial feature that precisely centres a limited volume insert in the vial neck. Available with graduated marking spots which can be written on with a #2 pencil.

Agilent	Velocity	Description
5183-4507	32011S-1232Z	2.0mL Clear Snap Top Vial 100 Pack - Deactivated
5183-4508	32011SE-1232	2.0mL Clear Wide Mouth Snap Top Vial with White Graduation Spot 100 Pack
5183-4509	32011SE-1232A	2.0mL Amber Wide Mouth Top Vial with White Graduation Spot 100 Pack
5182-0541	5250-11B	11mm Blue Snap Cap, PTFE/Silicone Liner 100 Pack
5182-3458	5240-11B	11mm Blue Snap Cap, PTFE/Butyl Liner 100 Pack




30211S-1232



30511CP-1232, 30111CP-1232



11mm Aluminium



11mm Snap Caps

Vials	
Description	Part Number
<b>Standard Vials</b>	
2.0mL Clear Snap Top Vial 100 Pack	32011S-1232
2.0mL Clear Snap Top Vial 100 Pack - Deactivated	32011S-1232Z
2.0mL Amber Snap Top Vial 100 Pack	32011S-1232A
2.0mL Amber Snap Top Vial 100 Pack - Deactivated	32011S-1232AZ
2.0mL Clear Snap Top Vial with White Graduation Spot 100 Pack	32011SE-1232
2.0mL Amber Standard Snap Top Vial with White Graduation Spot 100 Pack	3211SE-1232A
<b>Interlocked Vials with Fused glass Inserts</b>	
100µL Clear Interlocked™ Vial/Insert 100 Pack	30111S-1232
300µL Clear Interlocked™ Vial/Insert 100 Pack	30211S-1232
300µL Amber Interlocked™ Vial/Insert 100 Pack	30211S-1232A
300µL Clear Interlocked™ Vial/Insert with Marking Spot 100 Pack	30211SM-1232
300µL Amber Interlocked™ Vial/Insert with Marking Spot 100 Pack	3021SM-1232A
<b>Plastic Snap Top Vials</b>	
100µL Clear Polypropylene Crimp/Snap Vial 100 Pack	30111CP-1232
500µL Clear Polypropylene Crimp/Snap Vial 100 Pack	30511CP-1232
750µL Clear Polypropylene Crimp/Snap Vial 100 Pack	30711CP-1232

Description	Part Number				
	Silver	Blue	Gold	Red	Green
11mm Seal, 10mil PTFE 100 Pack	5110-11	-	-	-	-
11mm Seal, Clear PTFE/Red Rubber 100 Pack	5140-11	5140-11B	5140-11Y	5140-11R	5140-11G
11mm Seal, PTFE/Silicone 100 Pack	5150-11	5150-11B	5150-11Y	5150-11R	5150-11G
11mm Seal, PTFE/Silicone/PTFE 100 Pack	5160-11	5160-11B	5160-11Y	5160-11R	5160-11G
11mm Seal, PTFE/Silicone w/Slit 100 Pack	5170-11	-	-	-	-

Description	Part Number				
	Clear	Blue	Yellow	Red	Green
11mm Snap, 10mil PTFE 100 Pack	5210-11	5210-11B	5210-11Y	5210-11R	5210-11G
11mm Snap, Clear PTFE/Red Rubber 100 Pack	5240-11	5240-11B	5240-11Y	5240-11R	5240-11G
11mm Snap, PTFE/Silicone 100 Pack	5250-11	5250-11B	5250-11Y	5250-11R	5250-11G
11mm Snap, PTFE/Silicone/PTFE 100 Pack	5260-11	5260-11B	5260-11Y	5260-11R	5260-11G
11mm Snap, PTFE/Silicone w/Slit 100 Pack	5270-11	-	-	-	-

## 2.0mL Wide Mouth Snap Top Vials, Inserts and Caps

### GC Snap Top Caps™

GC Snap Top Caps™ were developed to handle the most volatile samples, providing a liquid airtight seal to the vial. Simply snap on the cap and an evaporation-proof seal of a crimp cap is achieved.



GC Snap Caps					
Description	Part Number				
	Clear	Blue	Yellow	Red	Green
11mm GC Snap Cap, 10mil PTFE	5210-11PP	-	-	-	-
11mm GC Snap Cap, PTFE/Butyl Rubber	5240-11PP	-	-	-	-
11mm GC Snap Cap, PTFE/Silicone	5250-11PP	-	-	-	-
11mm GC Snap Cap, PTFE/Silicone/PTFE	5260-11PP	-	-	-	-



4005BS-625, 4025BS-529,  
401TS-529

Inserts	
Description	Part Number
<b>Conical Glass Inserts</b>	
50µL Glass Insert with Bottom Spring 100 Pack	4005BS-625
50µL Glass Insert with Bottom Spring 100 Pack - Deactivated	4005BS-625Z
100µL Glass Insert with Glass Flange 100 Pack	401TF-529
100µL Glass Insert with Top Spring 100 Pack	401TS-529
250µL Glass Big Mouth Insert, No Spring 100 Pack	4025-629
250µL Glass Big Mouth Insert, No Spring 100 Pack - Deactivated	4025-629Z
250µL Glass Big Mouth Insert with Bottom Spring 100 Pack	4025BS-629
250µL Glass Big Mouth Insert with Bottom Spring 100 Pack - Deactivated	4025BS-629Z
<b>Flat Bottom Glass Inserts</b>	
350µL Glass Flat Bottom Insert 100 Pack	4035-630
350µL Glass Flat Bottom Insert 100 Pack - Deactivated	4035-630Z
350µL Glass Flat Bottom Insert with ID Ring 100 Pack	4035-630ID
350µL Glass Flat Bottom Insert with ID Ring 100 Pack - Deactivated	4035-630IDZ
<b>Conical Polypropylene Inserts</b>	
100µL Polypropylene Insert with Top Spring 100 Pack	401PTS-530
250µL Polypropylene Big Mouth Insert, No Spring 100 Pack	4025P-631
250µL Polypropylene Big Mouth Insert with Bottom Spring 100 Pack	4025PBS-631



4035-630ID, 4035-630



4025P-631, 4025PBS-631



4025BS-629, 401TS-529