

Headspace Vials and Caps



Left: Flat Bottom, Screw Top Vial
Right: Round Bottom, Crimp Top Vial




320020-2375, 310020-2346



320018R-2375, 310018R-2346



BiMetal Seals. Coloured Outer surface is Aluminium making them easy to crimp but still magnetic due to inner tin surface.



Agilent	Velocity	Description
5182-0837	320020-2375	20mL Headspace Vial, 23x75mm, 20mm Crimp (Flat Bottom)
5182-0838	310020-2346	10mL Headspace Vial, 23x46mm, 20mm Crimp (Flat Bottom)
5183-4474	320020A-2375	20mL Headspace Vial, 23x75mm, 20mm Crimp (Rounded Bottom)
5183-4475	310020A-2346	10mL Headspace Vial, 23x46mm, 20mm Crimp (Rounded Bottom)
5183-4477	5150-20	20mm Silver Seal, PTFE/Silicone Lined
5183-4480	5141PR-20	20mm Silver Seal, PTFE/Gray Butyl Molded – Pressure Release
5183-4478	5143PR-20	20mm Silver Seal, PTFE/Silicone Pressure Release

Flat Bottom Crimp Top Vials — Use in Agilent, Perkin Elmer and Shimadzu	
Description	Part Number
6mL Clear Headspace Vial, 22x38mm (for Perkin-Elmer), 20mm Crimp 100 Pack	36020-2238
10mL Clear Headspace Vial, 23x46mm, 20mm Crimp 100 Pack	310020-2346
20mL Clear Headspace Vial, 23x75mm, 20mm Crimp 100 Pack	320020-2375
27mL Clear Headspace Vial, 30x60mm (for Shimadzu), 20mm Crimp 100 Pack	327020-3060

Flat Bottom Screw Top Vials — Use in Agilent, Perkin Elmer and Shimadzu	
Description	Part Number
10mL Clear Flat Bottom Headspace Vial, 23x46mm, 18mm Thread 100 Pack	310018F-2346
10mL Amber Flat Bottom Headspace Vial, 23x46mm, 18mm Thread 100 Pack	31018F-2346A
20mL Clear Flat Bottom Headspace Vial, 23x75mm, 18mm Thread 100 Pack	320018F-2375
20mL Amber Flat Bottom Headspace Vial, 23x75mm, 18mm Thread 100 Pack	32018F-2375A

Round Bottom Crimp Top Vials — Use in Perkin Elmer, L.E.A.P. and Varian	
Description	Part Number
10mL Headspace Vial, 23x46mm (for Leap), 20mm Crimp 100 Pack	310020A-2346
20mL Headspace Vial, 23x75mm (for Perkin Elmer), 20mm Crimp 100 Pack	320020A-2375

Round Bottom Screw Top Vials — Use in Perkin Elmer, L.E.A.P. and Varian	
Description	Part Number
10mL Clear Rounded Bottom Headspace Vial, 23x46mm, 18mm Thread 100 Pack	310018R-2346
10mL Amber Rounded Bottom Headspace Vial, 23x46mm, 18mm Thread 100 Pack	31018R-2346A
20mL Clear Rounded Bottom Headspace Vial, 23x75mm, 18mm Thread 100 Pack	320018R-2375
20mL Amber Rounded Bottom Headspace Vial, 23x75mm, 18mm Thread 100 Pack	32018R-2375A

20mm Crimp Seals	
Description	Part Number
20mm Silver Seal, 10mil PTFE 100 Pack	5110-20
20mm Silver Seal, 0.131" PTFE/Gray Butyl Rubber 100 Pack	5140-20
20mm Silver Seal, 0.125" PTFE/Silicone 100 Pack	5150-20
20mm Silver Seal, 0.131" PTFE/Molded Gray Butyl Pressure Release 100 Pack	5141PR-20
20mm Silver Seal, 0.100" PTFE/Silicone Pressure Release 100 Pack	5143PR-20
Magnetic Seals	
20mm Silver Tin Plate Seal, 5mm Hole, 0.131" PTFE/Butyl Rubber 100 Pack	5144TP5-20
20mm Silver Tin Plate Seal, 8mm Hole, 0.131" PTFE/Butyl Rubber 100 Pack	5144TP8-20
20mm Silver Tin Plate Seal, 5mm Hole, 0.100" PTFE/Silicone 100 Pack	5145TP5-20
20mm Silver Tin Plate Seal, 8mm Hole, 0.100" PTFE/Silicone 100 Pack	5145TP8-20
20mm Red BiMetal Seal, 0.118" PTFE/Silicone 100 Pack	5150BM18-20R
20mm Blue BiMetal Seal, 0.118" PTFE/Silicone 100 Pack	5150BM18-20B

18mm Screw Thread Seals	
Description	Part Number
18mm Silver Metal Closure Only 100 Pack	5310M-18
18mm Silver Metal Closure with 0.050" PTFE/Silicone Line 100 Pack	5350-18M
18mm Silver Metal Closure with Blue 0.060" PTFE/Silicone Liner 100 Pack	5356B-18M
18mm Silver Metal Closure with Trans Blue 0.125" PTFE/Silicone Liner 100 Pack	5351TB-18M

20mm Stoppers	
Description	Part Number
20mm Butyl rubber Stoppers 100 Pack	5001-20ST