

## Vial Racks, Crimpers and Accessories

### Hand Operated Vial Crimper

Use to crimp vials and bottles in one quick, easy motion.

Built-in adjustable stop allows the user to repeatedly achieve an optimum seal.



Hand Operated Vial Crimper	
Description	Part Number
8mm Hand Operated Crimper	9300-08
11mm Hand Operated Crimper	9300-11
13mm Hand Operated Crimper	9300-13
20mm Hand Operated Crimper	9300-20

### Plier Type Decapper

Plier type decappers are a convenient, low cost hand tool which allows crimped aluminium seals to be removed quickly and easily from vials and bottles.



Plier Type Decapper	
Description	Part Number
8mm Plier Type Decapper	<b>9310-08</b>
11mm Plier Type Decapper	<b>9310-11</b>
13mm Plier Type Decapper	<b>9310-13</b>
20mm Plier Type Decapper	<b>9310-20</b>

### Hand Operated Vial Decapper

Hand operated crimp type decappers are used to remove the crimp seal from vials and bottles in one smooth motion.

Sizes range from 8mm to 20mm.



Hand Operated Vial Decapper	
Description	Part Number
8mm Hand Operated Decapper	<b>9320-08</b>
11mm Hand Operated Decapper	<b>9320-11</b>
13mm Hand Operated Decapper	<b>9320-13</b>
20mm Hand Operated Decapper	<b>9320-20</b>

### The Multi-Plug Safety Sleeve

Use it as a tool to make plug and shell vial assembly easier and safer. The safety sleeve is "Multi" useful because it is designed to work with 1mL, 2mL or 4mL shell vials. One Safety Sleeve works with all three sizes.



The Multi-Plug Safety Sleeve	
Description	Part Number
Multi-Plug Safety Sleeve Clear	<b>9200-81215</b>
Multi-Plug Safety Sleeve Black	<b>9200-81215B</b>

### Vial Storage Boxes

Rugged polypropylene box has durable three-point hinges and clasp for secure closure, molded grid lines on the lid, imprinted coordinates on the bottom of base, and imprinted and raised coordinates on each well. Holds standard 12mm (D) x 32mm (H) vials. Autoclavable. Assorted colors: blue, green, purple, yellow, and orange. Measures 141 x 151 x 57 mm. Pack of 5.



Vial Storage Boxes	
Description	Part Number
2mL Vial Storage Box Assorted Colours. 5 Pack	<b>VSS-1065</b>

### Vial Snap Racks

Snap Racks™ are constructed of durable glass reinforced polypropylene, with a centre well to conveniently hold caps, seals, septa, etc. Racks snap together securely to double your rack capacity.

Vial Storage Boxes	
Description	Part Number
96 Position Snap Rack™ White Rack for 8mm Vials	<b>9700-08</b>
96 Position Snap Rack™ Blue Rack for 8mm Vials	<b>9700-08B</b>
50 Position Snap Rack™ White Rack for 12mm Vials	<b>9700-12</b>
50 Position Snap Rack™ Blue Rack for 12mm Vials	<b>9700-12B</b>

